













# TRULENE

## MONOFILAMENT POLYPROPYLENE SUTURE












## EYELESS NEEDLED SUTURES

NEEDLE	Profile & Description	LENGTH	SIZE				
			6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
 <p>12mm 12mm Double Armed</p>	3/8 Circle Round Bodied	70 cm		SN 827			
 <p>12mm</p>	3/8 Circle Reverse Cutting	45 cm	SN 878		SN 879		
		70 cm		SN 880			
 <p>13mm 13mm Double Armed</p>	3/8 Circle Round Bodied	70 cm	SN 829				
 <p>15 mm</p>	3/8 Circle Slim Blade Cutting	70 cm	SN 823	SN 826			
 <p>16 mm</p>	3/8 Circle Round Bodied	70 cm		SN 881			
 <p>16 mm 16 mm Double Armed</p>	3/8 Circle Round Bodied	70 cm		SN 882			
 <p>16 mm</p>	3/8 Circle Cutting	70 cm			SN 870		
 <p>16 mm</p>	1/2 Circle Round Bodied	70 cm			SN 849		
 <p>17 mm</p>	1/2 Circle Taper Point	70 cm			SN 899		
		90 cm			SN 857		
 <p>17 mm 17 mm Double Armed</p>	1/2 Circle Taper Point	70 cm			SN 835	SN 897	SN 898
		90 cm			SN 853	SN 836	SN 837
 <p>20 mm</p>	1/2 Circle Round Bodied	70 cm				SN 887	
 <p>22 mm</p>	3/8 Circle Cutting	70 cm				SN 018	

# TRULENE

## MONOFILAMENT POLYPROPYLENE SUTURE

## EYELESS NEEDLED SUTURES

NEEDLE	Profile & Description	LENGTH	SIZE				
			4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
 <p>25 mm</p>	1/2 Circle Round Bodied	70 cm		SN 838			
		90 cm		SN 825			
 <p>25 mm 25 mm</p> <p>Double Armed</p>	1/2 Circle Round Bodied	70 cm		SN 822			
		90 cm		SN 852			
 <p>25 mm</p>	1/2 Circle Taper Point	70 cm			SN 807	SN 805	
		90 cm			SN 806		
 <p>25 mm 25 mm</p> <p>Double Armed</p>	1/2 Circle Taper Point	70 cm			SN 808		
		90 cm	SN 867		SN 877		
 <p>26 mm 26 mm</p> <p>Double Armed</p>	1/2 Circle Round Bodied	90 cm				SN 824	
 <p>30 mm</p>	1/2 Circle Round Bodied	70 cm			SN 841	SN 842	
 <p>30 mm HEAVY</p>	1/2 Circle Round Bodied	70 cm				SN 830	SN 840
 <p>30 mm HEAVY</p>	1/2 Circle Taper Point	90 cm				SN 894	
		100 cm					SN 896
 <p>40 mm HEAVY</p>	1/2 Circle Round Bodied	70 cm				SN 843	
 <p>45 mm HEAVY</p>	1/2 Circle Reverse Cutting	100 cm				SN 883	
 <p>60 mm</p>	Straight Cutting	70 cm		SN 888			