










TRULON

MONOFILAMENT POLYAMIDE SUTURE

EYELESS NEEDLED SUTURES





NEEDLE	Profile & Description	LENGTH	SIZE			
			6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)
12 mm 	3/8 Circle Reverse Cutting	70 cm		SN 3317		
15 mm 	3/8 Circle Slim Blade Cutting	70 cm		SN 3726		
16 mm 	3/8 Circle Cutting	35 cm		SN 3315		
		70 cm			SN 3318	SN 3321
26 mm 	3/8 Circle Cutting	70 cm	SN 3331			
26 mm 	3/8 Circle Reverse Cutting	70 cm				SN 3328

NEEDLE	Profile & Description	LENGTH	SIZE				
			3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
38 mm 	1/2 Circle Reverse Cutting	75 cm		SN 3303	SN 3304	SN 3305	SN 3306
40 mm 	1/2 Circle Round Bodied	150 cm LOOP			SN 3340		
40 mm HEAVY 	1/2 Circle Round Bodied	100 cm			SN 3346	SN 3347	
45 mm 	3/8 Circle Reverse Cutting	70 cm		SN 3336			

TRULON

MONOFILAMENT POLYAMIDE SUTURE

EYELESS NEEDLED SUTURES

NEEDLE	Profile & Description	LENGTH	SIZE			
			3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)
 <p>45 mm HEAVY</p>	1/2 Circle Reverse Cutting	100 cm				SN 3338
 <p>50 mm 50 mm Double Armed</p>	Straight Round Bodied	35 cm	SN 3307			
 <p>50 mm HEAVY</p>	1/2 Circle Round Bodied	150 cm LOOP				SN 3348
 <p>60 mm</p>	Straight Cutting	70 cm	SN 3388	SN 3390	SN 3389	

PRE-CUT SUTURES

LENGTH	SIZE			
	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
2 x 76 cm	S 903	S 904	S 905	S 906