











TRUSILK

BRAIDED SILK SUTURE




EYELESS NEEDLED SUTURES

NEEDLE	Profile & Description	LENGTH	SIZE				
			3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
25 mm 	1/2 Circle Round Bodied	76 cm	SN 5070				
25 mm 	1/2 Circle Taper Point	76 cm		SN 5670			
26 mm 	3/8 Circle Reverse Cutting	76 cm	SN 5028				
30 mm 	3/8 Circle Round Bodied	76 cm			SN 5017		
30 mm 	1/2 Circle Round Bodied	76 cm		SN 5333	SN 5334		
35 mm 	1/2 Circle Round Bodied	76 cm	SN 5083				
38 mm 	1/2 Circle Reverse Cutting	75 cm		SN 5003	SN 5004	SN 5005	SN 5006
40 mm 	Straight Round Bodied	76 cm		SN 5256			
45 mm 	3/8 Circle Reverse Cutting	76 cm		SN 5036			
50 mm 	Straight Round Bodied	76 cm		SN 5208			

TRUSILK

BRAIDED SILK SUTURE

EYELESS NEEDLED SUTURES

NEEDLE	Profile & Description	LENGTH	SIZE		
			2-0 (3)	1-0 (3.5)	1 (4)
 <p>50 mm</p>	1/2 Circle Reverse Cutting	76 cm	SN 5052		
 <p>60 mm</p>	Straight Cutting	76 cm	SN 5290		
 <p>60 mm</p>	3/8 Circle Cutting	76 cm			SN 5062

PRE - CUT SUTURES

LENGTH	PACKAGING	SIZE							
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
2 X 75 cm	STERILE FOIL OVERWRAP PACKS	S 209	S 210	S 211	S 212	S 213	S 214	S 215	S 216

REELS

LENGTH	SIZE							
	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	1-0 (3.5)	1 (4)	2 (5)
25 METERS A Box of 6 Reels	B 819	B 820	B 821	B 822	B 823	B 824	B 825	B 826